

AHNK Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP06118
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human protein . at AA range: 5720-5800
Mol wt	N/A
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	IHC-p, IF
Concentration	1 mg/ml
Full name	Neuroblast differentiation-associated protein AHNAK
Synonyms	Neuroblast differentiation-associated protein AHNAK ;Desmoyokin

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

May be required for neuronal cell differentiation.,Contains 1 PDZ (DHR) domain.,subunit:Interacts with DYSF; the interaction is direct and Ca(2+)-independent.,

Recommended Dilution

IHC-p: 1: 50 - 1: 300

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com