

ALKB4 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP14186
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Synthesized peptide derived from part region of human protein
Mol wt	N/A
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Probable alpha-ketoglutarate-dependent dioxygenase ABH4
Synonyms	Probable alpha-ketoglutarate-dependent dioxygenase ABH4 ;EC 1.14.11.-.;Alkylated DNA repair protein alkB homolog 4

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

sequence Wrong choice of CDS.,Belongs to the alkB family.,


Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images



Western blot analysis of lysates from SW480 cells, primary antibody was diluted at 1:1000, 4°C overnight

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com