

P22 Monoclonal Antibody

Description

Product type	Antibody
Code	BT-MCA3323
Host	Mouse
Isotype	Mouse IgG1
Size	100µL, 50µL
Immunogen	Purified recombinant fragment of human P22 (AA: 29-177) expressed in E. Coli.
Mol wt	N/A
Species reactivity	Others
Clonality	Monoclonal
Recommended application	Others
Concentration	N/A
Full name	N/A
Synonyms	KP177R

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

N/A

Recommended Dilution

Not yet tested in other applications.

Images

No images.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term.