

DNAL4 Monoclonal Antibody

Description

Product type	Antibody
Code	BT-MCA3907
Host	Mouse
Isotype	Mouse IgG1
Size	100µL, 50µL
Immunogen	Purified recombinant fragment of human DNAL4 (AA: 1-105) expressed in E. Coli.
Mol wt	12kDa
Species reactivity	Human
Clonality	Monoclonal
Recommended application	IHC
Concentration	N/A
Full name	N/A
Synonyms	PIG27

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

DNAL4 is a component of the dynein motor complex

Recommended Dilution

WB: 1:500 - 1:2000

IHC-p: 1:200 - 1:1000

ELISA: 1:10000

Not yet tested in other applications.

Images

No images.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term.

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com