

Influenza A virus Nucleoprotein Monoclonal Antibody

Description

Product type	Antibody
Code	BT-MCA3986
Host	Mouse
Isotype	Mouse IgG1
Size	100μL, 50μL
Immunogen	Purified recombinant fragment of Influenza A virus nucleoprotein (AA: 2-498) expressed in E. Coli strain BL21(DE3)
Mol wt	N/A
Species reactivity	Others
Clonality	Monoclonal
Recommended application	N/A
Concentration	N/A
Full name	N/A
Synonyms	N/A

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

N/A


Recommended Dilution

WB: 1:500 - 1:2000

ELISA: 1:10000

Not yet tested in other applications.

Images



Western blot analysis using Influenza A virus Nucleoprotein mouse mAb against full-length recombinant Influenza A virus Nucleoprotein.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term.