

ALK Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP00370
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	The antiserum was produced against synthesized peptide derived from human ALK. AA range:1473-1522
Mol wt	176428
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	ALK Antibody
Synonyms	ALK; ALK tyrosine kinase receptor; Anaplastic lymphoma kinase; CD antigen CD246

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

ALK encodes a receptor tyrosine kinase, which belongs to the insulin receptor superfamily. ALK tyrosine kinase receptor comprises an extracellular domain, an hydrophobic stretch corresponding to a single pass transmembrane region, and an intracellular kinase domain. It plays an important role in the development of the brain and exerts its effects on specific neurons in the nervous system. This gene has been found to be rearranged, mutated, or amplified in a series of tumours including anaplastic large cell lymphomas, neuroblastoma, and non-small cell lung cancer. The chromosomal rearrangements are the most common genetic alterations in this gene, which result in creation of multiple fusion genes in tumourigenesis, including ALK (chromosome 2)/EML4 (chromosome 2), ALK/RANBP2 (chromosome 2), ALK/ATIC (chromosome 2), ALK/TFG (chromosome 3), ALK/NPM1 (chromosome 5), ALK/SQSTM1 (chromosome 5), ALK/KIF5B (chromosome 10), ALK/CLTC (chromosome 17), ALK/TPM4 (chromosome 19), and ALK/MSN (chromosome X).


Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 5000

Not yet tested in other applications.

Images



Western blot analysis of lysates from LOVO cells, using ALK Antibody. The lane on the right is blocked with the synthesized peptide.

Storage

-20°C for one year