

HHEX Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP09883
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Synthesized peptide derived from part region of human protein
Mol wt	N/A
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Hematopoietically-expressed homeobox protein HHEX
Synonyms	Hematopoietically-expressed homeobox protein HHEX ;Homeobox protein HEX;Homeobox protein PRH

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

This gene encodes a member of the homeobox family of transcription factors, many of which are involved in developmental processes. Expression in specific hematopoietic lineages suggests that this protein may play a role in hematopoietic differentiation.


Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images



Western blot analysis of lysates from A431 cells, primary antibody was diluted at 1:1000, 4°C overnight

Storage

-20°C for 1 year