

POK7 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP12054
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human POK7
Mol wt	89320
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	HERV-K_1q23.3 provirus ancestral Pol protein
Synonyms	HERV-K_1q23.3 provirus ancestral Pol protein ;HERV-K;C1a Pol protein;HERV-K110 Pol protein;HERV-K18 Pol protein [Includes: Reverse transcriptase ;EC 2.7.7.49; Ribonuclease H ;RNase H;EC 3.1.2

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com