

NEUA Polyclonal Antibody

Description

Product type	Antibody
Code	BT-AP15813
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human C-terminal NEUA
Mol wt	48kD
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	IHC, WB
Concentration	1 mg/ml
Full name	N/A
Synonyms	N/A

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

N/A


Recommended Dilution

WB: 1:500-2000

IHC-p: 1:50-200

Not yet tested in other applications.

Images



Western Blot analysis of HeLa cell lysis, using primary antibody at 1:1000 dilution. Secondary antibody(AS0415) was diluted at 1:10000

Storage

-20°C for one year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com