

DNAJB9 Polyclonal Antibody

Description

Product type	Antibody
Code	BT-AP15816
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human DNAJB9
Mol wt	25kD
Species reactivity	Human,Mouse,Rat
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	N/A
Synonyms	DnaJ homolog subfamily B member 9 (Microvascular endothelial differentiation gene 1 protein) (Mdg-1)

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

N/A

Recommended Dilution

WB: 1:500-2000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for one year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com