

CD334 Monoclonal Antibody

Description

| | |
|--------------------------------|---|
| Product type | Antibody |
| Code | BT-MCA2154 |
| Host | Mouse |
| Isotype | Mouse IgG1 |
| Size | 100µL, 50µL |
| Immunogen | Purified recombinant fragment of human CD334 (AA: extra 22-369) expressed in E. Coli. |
| Mol wt | 88kDa |
| Species reactivity | Others |
| Clonality | Monoclonal |
| Recommended application | WB,FCM |
| Concentration | N/A |
| Full name | N/A |
| Synonyms | FGFR4;TKF;JTK2 |

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

The protein encoded by this gene is a tyrosine kinase and cell surface receptor for fibroblast growth factors. The encoded protein is involved in the regulation of several pathways, including cell proliferation, cell differentiation, cell migration, lipid metabolism, bile acid biosynthesis, vitamin D metabolism, glucose uptake, and phosphate homeostasis. This protein consists of an extracellular region, composed of three immunoglobulin-like domains, a single hydrophobic membrane-spanning segment, and a cytoplasmic tyrosine kinase domain. The extracellular portion interacts with fibroblast growth factors, setting in motion a cascade of downstream signals, ultimately influencing mitogenesis and differentiation. [provided by RefSeq, Aug 2017]

Recommended Dilution


WB: 1:500 - 1:2000

FCM: 1:200 - 1:400


ELISA: 1:10000

Not yet tested in other applications.


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)




Western blot analysis using CD334 mAb against human CD334 (AA: extra 22-369) recombinant protein. (Expected MW is 64.5 kDa)



Western blot analysis using CD334 mAb against HEK293 (1) and CD334 (AA: extra 22-369)-hlgGFc transfected HEK293 (2) cell lysate.



Western blot analysis using CD334 mouse mAb against K562 (1), MCF-7 (2), COS7 (3), and PC-3 (4) cell lysate.



Flow cytometric analysis of HL-60 cells using CD334 mouse mAb (green) and negative control (red).

Storage

Store at 4°C short term. Aliquot and store at -20°C long term.

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com