

KRT2 Monoclonal Antibody

Description

Product type Antibody

Code BT-MCA2205

Host Mouse

 Isotype
 Mouse IgG2a

 Size
 100μL, 50μL

Immunogen Purified recombinant fragment of human KRT2 (AA: 177-490) expressed in E. Coli.

Mol wt 65.4kDa

Species reactivity Others

Clonality Monoclonal

Recommended application IHC,ICC

Concentration N/A
Full name N/A

Synonyms K2e;KRTE;CK-2e;KRT2A;KRT2E

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background


The protein encoded by this gene is a member of the keratin gene family. The type II cytokeratins consist of basic or neutral proteins which are arranged in pairs of heterotypic keratin chains coexpressed during differentiation of simple and stratified epithelial tissues. This type II cytokeratin is expressed largely in the upper spinous layer of epidermal keratinocytes and mutations in this gene have been associated with bullous congenital ichthyosiform erythroderma. The type II cytokeratins are clustered in a region of chromosome 12q12-q13.

Recommended Dilution

WB: 1:500 - 1:2000 IHC-p: 1:200 - 1:1000 ICC: 1:200 - 1:1000 ELISA: 1:10000

Not yet tested in other applications.


Images




26-

Black line: Control Antigen (100 ng); Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)


Western blot analysis using KRT2 mAb against human KRT2 (AA: 177-490) recombinant protein. (Expected MW is 52.4 kDa)




Western blot analysis using KRT2 mAb against HEK293-6e (1) and KRT2 (AA: 177-490)-hIgGFc transfected HEK293-6e (2) cell lysate.




Western blot analysis using KRT2 mouse mAb against A431 (1) cell lysate.



Immunofluorescence analysis of Hela cells using KRT2 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)



Immunohistochemical analysis of paraffin-embedded ovarian cancer tissues using KRT2 mouse mAb with DAB staining.



 $Immunohistochemical\ analysis\ of\ paraffin-embedded\ breast\ cancer\ tissues\ using\ KRT2\ mouse\ mAb$ with DAB staining.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term.

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com