

ERCC1 Monoclonal Antibody

Description

Product type	Antibody
Code	BT-MCA2224
Host	Mouse
Isotype	Mouse IgG1
Size	100μL, 50μL
Immunogen	Purified recombinant fragment of human ERCC1 (AA: 1-297) expressed in E. Coli.
Mol wt	32.5kDa
Species reactivity	Others
Clonality	Monoclonal
Recommended application	WB,IHC,FCM
Concentration	N/A
Full name	N/A
Synonyms	UV20;COFS4;RAD10

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

The product of this gene functions in the nucleotide excision repair pathway, and is required for the repair of DNA lesions such as those induced by UV light or formed by electrophilic compounds including cisplatin. The encoded protein forms a heterodimer with the XPF endonuclease (also known as ERCC4), and the heterodimeric endonuclease catalyzes the 5' incision in the process of excising the DNA lesion. The heterodimeric endonuclease is also involved in recombinational DNA repair and in the repair of inter-strand crosslinks. Mutations in this gene result in cerebrooculofacioskeletal syndrome, and polymorphisms that alter expression of this gene may play a role in carcinogenesis. Multiple transcript variants encoding different isoforms have been found for this gene. The last exon of this gene overlaps with the CD3e molecule, epsilon associated protein gene on the opposite strand.

Recommended Dilution

WB: 1:500 - 1:2000


IHC-p: 1:200-1:1000

FCM: 1:200-1:400


ELISA: 1:10000

Not yet tested in other applications.


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)




Western blot analysis using ERCC1 mAb against human ERCC1 (AA: 1-297) recombinant protein. (Expected MW is 35.5 kDa)




Western blot analysis using ERCC1 mAb against HEK293-6e (1) and ERCC1 (AA: 1-297)-hIgGFc transfected HEK293 (2) cell lysate.




Western blot analysis using ERCC1 mouse mAb against NIH/3T3 (1), MCF-7 (2), HeLa (3), SK-Br-3 (4), HepG2 (5), Raji (6), PC-3 (7), and K562 (8) cell lysate.



Flow cytometric analysis of HeLa cells using ERCC1 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded lung cancer tissues using ERCC1 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded kidney tissues using ERCC1 mouse mAb with DAB staining.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term.

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com