

ELK1 Monoclonal Antibody

Description

Product type Antibody

Code BT-MCA2234

Host Mouse

 Isotype
 Mouse IgG1

 Size
 100µL, 50µL

Immunogen Purified recombinant fragment of ELK1 expressed in E. Coli.

Mol wt N/A

Species reactivity Human

Clonality Monoclonal

Recommended application Others

Concentration N/A

Full name N/A

Synonyms ELK1

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

N/A

Recommended Dilution

WB: 1:500 - 1:2000 ELISA: 1:10000

Not yet tested in other applications.

Images



 $We stern \ blot \ analysis \ using \ ELK1 \ mouse \ mAb \ against \ truncated \ ELK1 \ recombinant \ protein \ (1) \ and \ K562 \ cell \ ly sate \ (2).$

Storage

Store at 4°C short term. Aliquot and store at -20°C long term.