

MESP2 Monoclonal Antibody

Description

Product type Antibody

Code BT-MCA2431

Host Mouse

 Isotype
 Mouse IgG1

 Size
 100μL, 50μL

Immunogen Purified recombinant fragment of human MESP2 (AA: 37-94) expressed in E. Coli.

Mol wt 41.8kDa

Species reactivity Human

Clonality Monoclonal

Recommended application Others

Concentration N/A

Full name N/A

Synonyms SCDO2;bHLHc6

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

N/A

Recommended Dilution

WB: 1:500 - 1:2000 ELISA: 1:10000


Not yet tested in other applications.

Images




 $We stern \ blot \ analysis \ using \ MESP2 \ mAb \ against \ human \ MESP2 \ recombinant \ protein. \ (Expected \ MW)$

is 31.4 kDa)



Western blot analysis using MESP2 mAb against HEK293 (1) and MESP2 (AA: 37-94)-hIgGFc transfected HEK293 (2) cell lysate.



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);

Storage

Store at 4°C short term. Aliquot and store at -20°C long term.

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com