

KDM1A Monoclonal Antibody

Description

Product type	Antibody
Code	BT-MCA2603
Host	Mouse
Isotype	Mouse IgG1
Size	100μL, 50μL
Immunogen	Purified recombinant fragment of human KDM1A (AA: 55-263) expressed in E. Coli.
Mol wt	93kDa
Species reactivity	Human,Monkey
Clonality	Monoclonal
Recommended application	WB,IHC,ICC,FCM
Concentration	N/A
Full name	N/A
Synonyms	AOF2;CPRF;KDM1;LSD1;BHC110
This product is for proceeding only not for use in human thempositie or discretis procedure	

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.


Background

N/A

Recommended Dilution


WB: 1:500 - 1:2000 IHC-p: 1:200 - 1:1000 ICC: 1:200 - 1:1000 FCM: 1:200 - 1:400 ELISA: 1:10000 Not yet tested in other applications.

Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)

Western blot analysis using KDM1A mAb against human KDM1A (AA: 55-263) recombinant protein. (Expected MW is 60.3 kDa)



Storage

Store at 4°C short term. Aliquot and store at -20°C long term.

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com