

OTUD4 Monoclonal Antibody

Description

Product type	Antibody
Code	BT-MCA2768
Host	Mouse
Isotype	Mouse IgG1
Size	100 μ L, 50 μ L
Immunogen	Purified recombinant fragment of human OTUD4 (AA: 815-1049) expressed in E. Coli.
Mol wt	124kDa
Species reactivity	Human,Rat
Clonality	Monoclonal
Recommended application	WB,FCM
Concentration	N/A
Full name	N/A
Synonyms	HIN1;DUBA6;HSHIN1

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

N/A

Recommended Dilution


WB: 1:500 - 1:2000

FCM: 1:200 - 1:400

ELISA: 1:10000

Not yet tested in other applications.

Images



Flow cytometric analysis of HeLa cells using OTUD4 mouse mAb (green) and negative control (red).

Storage

Store at 4°C short term. Aliquot and store at -20°C long term.

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com