

HIST2H4A(20Me3) Monoclonal Antibody

Description

Product type	Antibody
Code	BT-MCA2941
Host	Mouse
Isotype	Mouse IgG1
Size	100μL, 50μL
Immunogen	Synthesized peptide of human HIST2H4A (AA: GGAKRHRK(Me3)VLRDNIQ) .
Mol wt	11.4kDa
Species reactivity	Human
Clonality	Monoclonal
Recommended application	FCM
Concentration	N/A
Full name	N/A
Synonyms	H4;H4/n;H4F2;H4FN;FO108;HIST2H4

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

N/A


Recommended Dilution

FCM: 1:200 - 1:400


ELISA: 1:10000

Not yet tested in other applications.

Images



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Flow cytometric analysis of Raji cells using HIST2H4A(20Me3) mouse mAb (green) and negative control (red).

Storage

Store at 4°C short term. Aliquot and store at -20°C long term.

