

# CD283 Monoclonal Antibody

# Description

Product type Antibody

Code BT-MCA2970

Host Mouse

 Isotype
 Mouse IgG1

 Size
 100μL, 50μL

Immunogen Purified recombinant fragment of human CD283 (AA: 24-217) expressed in E. Coli.

Mol wt 103.8kDa

Species reactivity Human

Clonality Monoclonal

Recommended application FCM

Concentration N/A

Full name N/A

Synonyms TLR3 ;IIAE2

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

# Background


The protein encoded by this gene is a member of the Toll-like receptor (TLR) family which plays a fundamental role in pathogen recognition and activation of innate immunity. TLRs are highly conserved from Drosophila to humans and share structural and functional similarities. They recognize pathogen-associated molecular patterns (PAMPs) that are expressed on infectious agents, and mediate the production of cytokines necessary for the development of effective immunity. The various TLRs exhibit different patterns of expression. This receptor is most abundantly expressed in placenta and pancreas, and is restricted to the dendritic subpopulation of the leukocytes. It recognizes dsRNA associated with viral infection, and induces the activation of NF-kappaB and the production of type I interferons. It may thus play a role in host defense against viruses. Use of alternative polyadenylation sites to generate different length transcripts has been noted for this gene. [provided by RefSeq, Jul 2008]

#### Recommended Dilution


WB: 1:500 - 1:2000 FCM: 1:200 - 1:400 ELISA: 1:10000

Not yet tested in other applications.


# **Images**




Black line: Control Antigen (100 ng); Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)



Western blot analysis using CD283 mAb against human CD283 (AA: 24-217) recombinant protein. (Expected MW is 47.5 kDa)



Western blot analysis using CD283 mAb against HEK293 (1) and CD283 (AA: 24-217)-hIgGFc transfected HEK293 (2) cell lysate.



Flow cytometric analysis of HL-60 cells using CD283 mouse mAb (green) and negative control (red).

# Storage

Store at 4°C short term. Aliquot and store at -20°C long term.

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com