

CD3D Monoclonal Antibody

Description

Product type	Antibody
Code	BT-MCA3037
Host	Mouse
Isotype	Mouse IgG2b
Size	100µL, 50µL
Immunogen	Purified recombinant fragment of human CD3D (AA: 127-171) expressed in E. Coli.
Mol wt	18.9kDa
Species reactivity	Others
Clonality	Monoclonal
Recommended application	WB,FCM
Concentration	N/A
Full name	N/A
Synonyms	T3D;IMD19;CD3-DELTA

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

The protein encoded by this gene is part of the T-cell receptor/CD3 complex (TCR/CD3 complex) and is involved in T-cell development and signal transduction. The encoded membrane protein represents the delta subunit of the CD3 complex, and along with four other CD3 subunits, binds either TCR alpha/beta or TCR gamma/delta to form the TCR/CD3 complex on the surface of T-cells. Defects in this gene are a cause of severe combined immunodeficiency autosomal recessive T-cell-negative/B-cell-positive/NK-cell-positive (SCIDBNK). Two transcript variants encoding different isoforms have been found for this gene. Other variants may also exist, but the full-length nature of their transcripts has yet to be defined. [provided by RefSeq, Feb 2009]

Recommended Dilution


WB: 1:500 - 1:2000

FCM: 1:200-1:400


ELISA: 1:10000

Not yet tested in other applications.


Images




Black line: Control Antigen (100 ng); Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)




Western blot analysis using CD3D mAb against human CD3D (AA: 127-171) recombinant protein. (Expected MW is 18.6 kDa)




Western blot analysis using CD3D mAb against HEK293-6e (1) and CD3D (AA: 127-171)-hIgGFc transfected HEK293-6e (2) cell lysate.



Western blot analysis using CD3D mouse mAb against HL-60 (1), and U937 (2) cell lysate.



Flow cytometric analysis of Jurkat cells using CD3D mouse mAb (green) and negative control (red).



Flow cytometric analysis of THP-1 cells using CD3D mouse mAb (green) and negative control (red).

Storage

Store at 4°C short term. Aliquot and store at -20°C long term.

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com