

E7 Monoclonal Antibody

Description

| | |
|--------------------------------|---|
| Product type | Antibody |
| Code | BT-MCA3072 |
| Host | Mouse |
| Isotype | Mouse IgG1 |
| Size | 100 μ L, 50 μ L |
| Immunogen | Purified recombinant fragment of E7 expressed in E. Coli. |
| Mol wt | N/A |
| Species reactivity | Others |
| Clonality | Monoclonal |
| Recommended application | Others |
| Concentration | N/A |
| Full name | N/A |
| Synonyms | E7, HpV16gp2 |

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

E7: E7 transforming protein (Human papillomavirus type 16).


Recommended Dilution

WB: 1:500 - 1:2000

ELISA: 1:10000

Not yet tested in other applications.

Images



Western blot analysis using E7 mouse mAb against truncated GST-E7 recombinant protein (1).

Storage

Store at 4°C short term. Aliquot and store at -20°C long term.

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com