

D117L Monoclonal Antibody

Description

Product type	Antibody
Code	BT-MCA3093
Host	Mouse
Isotype	Mouse IgG1
Size	100 μ L, 50 μ L
Immunogen	Purified recombinant fragment of human D117L (AA: 1-117) expressed in E. Coli.
Mol wt	N/A
Species reactivity	Others
Clonality	Monoclonal
Recommended application	WB
Concentration	N/A
Full name	N/A
Synonyms	p17

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background


structural protein p17 [African swine fever virus]

Recommended Dilution

WB: 1:500 - 1:2000

Not yet tested in other applications.

Images



Western blot analysis using D117L mAb against human D117L (AA: 1-117) recombinant protein. lane 1 :(100 ng); lane 2 :(50 ng); lane 3 :(25 ng); lane 4 :(10 ng); lane 5 :(2.5 ng); (Expected MW is 40 kDa)

Storage

Store at 4°C short term. Aliquot and store at -20°C long term.

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com