

SARS-mpm Monoclonal Antibody

Description

Product type	Antibody
Code	BT-MCA3150
Host	Mouse
Isotype	Mouse IgG2b
Size	100μL, 50μL
Immunogen	Purified recombinant fragment of SARS-mpm expressed in E. Coli.
Mol wt	N/A
Species reactivity	Others
Clonality	Monoclonal
Recommended application	Others
Concentration	N/A
Full name	N/A
Synonyms	N/A

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

N/A


Recommended Dilution

WB: 1:500 - 1:2000


ELISA: 1:10000

Not yet tested in other applications.

Images



Black line: Control Antigen (100 ng); Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)



Western blot analysis using SARS-mpm mAb against SARS-mpm recombinant protein. lane 1: 100 ng lane 2: 50 ng lane 3: 25 ng lane 4: 5 ng lane 5: 2.5 ng (Expected MW is 34.9 kDa)

Storage

Store at 4°C short term. Aliquot and store at -20°C long term.