

OVA Monoclonal Antibody

Description

Product type	Antibody
Code	BT-MCA3236
Host	Mouse
Isotype	Mouse IgG1
Size	100μL, 50μL
Immunogen	Ovalbumin
Mol wt	N/A
Species reactivity	Others
Clonality	Monoclonal
Recommended application	WB
Concentration	N/A
Full name	N/A
Synonyms	N/A

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Ovalbumin (OVA) is used in variety of laboratories products which contain OVA as a carrier protein. To raise antibodies against small non-immunogenic molecules, called haptens, it is often necessary to couple them to carrier protein, such as OVA. This antibody is a good reagent that may be used for the development of detection assay in ELISA or western blot to test the remaining OVA for clinical products in Vivo or research products In Vitro.


Recommended Dilution

WB: 1:500 - 1:2000

ELISA: 1:10000

Not yet tested in other applications.

Images



Western blot analysis using OVA mouse mAb against OVA protein.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term.