

LRP6 Monoclonal Antibody

Description

Product type Antibody

Code BT-MCA3360

Host Mouse

 Isotype
 Mouse IgG1

 Size
 100µL, 50µL

Immunogen Purified recombinant fragment of human LRP6 (AA: extra 1198-1370) expressed in E. Coli.

Mol wt 180kDa

Species reactivity Human

Clonality Monoclonal

Recommended application FCM

Concentration N/A
Full name N/A

Synonyms ADCAD2;STHAG7

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background


This gene encodes a member of the low density lipoprotein (LDL) receptor gene family. LDL receptors are transmembrane cell surface proteins involved in receptor-mediated endocytosis of lipoprotein and protein ligands. The protein encoded by this gene functions as a receptor or, with Frizzled, a co-receptor for Wnt and thereby transmits the canonical Wnt/beta-catenin signaling cascade. Through its interaction with the Wnt/beta-catenin signaling cascade this gene plays a role in the regulation of cell differentiation, proliferation, and migration and the development of many cancer types. This protein undergoes gamma-secretase dependent RIP- (regulated intramembrane proteolysis) processing but the precise locations of the cleavage sites have not been determined.

Recommended Dilution

WB: 1:500 - 1:2000 FCM: 1:200 - 1:400 ELISA: 1:10000


Not yet tested in other applications.

Images




Black line: Control Antigen (100 ng); Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)

Western blot analysis using LRP6 mAb against human LRP6 (AA: extra 1198-1370) recombinant protein. (Expected MW is 45 kDa)



Western blot analysis using LRP6 mAb against HEK293 (1) and LRP6 (AA: extra 1198-1370)-hIgGFc transfected HEK293 (2) cell lysate.



Flow cytometric analysis of HL-60 cells using LRP6 mouse mAb (green) and negative control (red).

Storage

Store at 4°C short term. Aliquot and store at -20°C long term.

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com