

KDM3A Monoclonal Antibody

Description

Product type	Antibody
Code	BT-MCA3449
Host	Mouse
Isotype	Mouse IgG1
Size	100µL, 50µL
Immunogen	Purified recombinant fragment of human KDM3A (AA: 311-445) expressed in E. Coli.
Mol wt	147kDa
Species reactivity	Human
Clonality	Monoclonal
Recommended application	FCM
Concentration	N/A
Full name	N/A
Synonyms	TSGA;JMJD1;JHDM2A;JHMD2A;JMJD1A

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

N/A

Recommended Dilution


WB: 1:500 - 1:2000

FCM: 1:200 - 1:400


ELISA: 1:10000

Not yet tested in other applications.


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using KDM3A mAb against human KDM3A (AA: 311-445) recombinant protein. (Expected MW is 40.3 kDa)



Western blot analysis using KDM3A mAb against HEK293 (1) and KDM3A (AA: 311-445)-hIgGFc transfected HEK293 (2) cell lysate.



Flow cytometric analysis of HeLa cells using KDM3A mouse mAb (green) and negative control (red).

Storage

Store at 4°C short term. Aliquot and store at -20°C long term.

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com