

MAPK10 Monoclonal Antibody

Description

Product type	Antibody
Code	BT-MCA3464
Host	Mouse
Isotype	Mouse IgG2a
Size	100µL, 50µL
Immunogen	Purified recombinant fragment of human MAPK10 (AA: 180-329) expressed in E. Coli.
Mol wt	52.6kDa
Species reactivity	Human,Rat
Clonality	Monoclonal
Recommended application	WB
Concentration	N/A
Full name	N/A
Synonyms	JNK3;JNK3A;PRKM10;SAPK1b;p493F12;p54bSAPK

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

N/A


Recommended Dilution

WB: 1:500 - 1:2000


ELISA: 1:10000

Not yet tested in other applications.


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using MAPK10 mAb against human MAPK10 (AA: 180-329) recombinant protein. (Expected MW is 43.1 kDa)



Western blot analysis using MAPK10 mAb against HEK293 (1) and MAPK10 (AA: 180-329)-hIgGfC transfected HEK293 (2) cell lysate.



Western blot analysis using MAPK10 mouse mAb against HEK293 (1), HeLa (2), SK-N-SH (3), MCF-7 (4), Jurkat (5), and C6 (6) cell lysate.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term.

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com