

CD281 Monoclonal Antibody

Description

Product type Antibody

Code BT-MCA3593

Host Mouse

 Isotype
 Mouse IgG1

 Size
 100μL, 50μL

Immunogen Purified recombinant fragment of human CD281 (AA: extra 172-425) expressed in E. Coli.

Mol wt 90.2kDa

Species reactivity Others

Clonality Monoclonal

Recommended application IHC,ICC,FCM

Concentration N/A
Full name N/A

Synonyms TLR1;TIL;rsc786;TIL. LPRS5

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background


The protein encoded by this gene is a member of the Toll-like receptor (TLR) family which plays a fundamental role in pathogen recognition and activation of innate immunity. TLRs are highly conserved from Drosophila to humans and share structural and functional similarities. They recognize pathogen-associated molecular patterns (PAMPs) that are expressed on infectious agents, and mediate the production of cytokines necessary for the development of effective immunity. The various TLRs exhibit different patterns of expression. This gene is ubiquitously expressed, and at higher levels than other TLR genes. Different length transcripts presumably resulting from use of alternative polyadenylation site, and/or from alternative splicing, have been noted for this gene.

Recommended Dilution

WB: 1:500 - 1:2000 IHC-p: 1:200 - 1:1000 ICC: 1:200 - 1:1000 FCM: 1:200 - 1:400 ELISA: 1:10000


Not yet tested in other applications.

Images




34-26Black line: Control Antigen (100 ng); Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)


Western blot analysis using CD281 mAb against human CD281 (AA: extra 172-425) recombinant protein. (Expected MW is 55.1 kDa)




Western blot analysis using CD281 mAb against HEK293-6e (1) and CD281 (AA: extra 172-425)-hIgGFc transfected HEK293-6e (2) cell lysate.




Immunofluorescence analysis of Hela cells using CD281 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)



Flow cytometric analysis of THP-1 cells using CD281 mouse mAb (green) and negative control (red).



 $Immun ohistochemical \ analysis \ of \ paraffin-embedded \ lung \ tissues \ using \ CD281 \ mouse \ mAb \ with \ DAB \ staining.$



Immunohistochemical analysis of paraffin-embedded rectum cancer tissues using CD281 mouse mAb with DAB staining.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term.

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com