

## PKD1 Monoclonal Antibody

### Description

<b>Product type</b>	Antibody
<b>Code</b>	BT-MCA3644
<b>Host</b>	Rat
<b>Isotype</b>	Rat IgG1
<b>Size</b>	100μL, 50μL
<b>Immunogen</b>	Purified recombinant fragment of mouse PKD1 (AA: 3682-3873) expressed in E. Coli.
<b>Mol wt</b>	466kDa
<b>Species reactivity</b>	Human
<b>Clonality</b>	Monoclonal
<b>Recommended application</b>	FCM
<b>Concentration</b>	N/A
<b>Full name</b>	N/A
<b>Synonyms</b>	PC1;Mflj00285

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

### Background

N/A


### Recommended Dilution

FCM: 1:200 - 1:400


ELISA: 1:10000

Not yet tested in other applications.

### Images



Flow cytometric analysis of HeLa cells using PKD1 Rat mAb (green) and negative control (red).



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);

### Storage

Store at 4°C short term. Aliquot and store at -20°C long term.

