

ALDH2 Monoclonal Antibody

Description

Product type Antibody

Code BT-MCA3686

Host Mouse

 Isotype
 Mouse IgG1

 Size
 100μL, 50μL

Immunogen Purified recombinant fragment of human ALDH2 (AA: 317-517) expressed in E. Coli.

Mol wt 56.3kDa

Species reactivity Human, Mouse, Rat

Clonality Monoclonal

Recommended application WB,IHC,ICC,FCM

Concentration N/A
Full name N/A

Synonyms ALDM;ALDHI;ALDH-E2

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

N/A

Recommended Dilution


WB: 1:500 - 1:2000 IHC-p: 1:200 - 1:1000

ICC: 1:150

FCM: 1:200 - 1:400 ELISA: 1:10000

Not yet tested in other applications.

Images



Western blot analysis using ALDH2 mouse mAb against HepG2 (1), A549 (2) cell lysate and Rat liver (3), Mouse liver (4), Mouse brain (5) tissue lysate.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term.