

## D117L Monoclonal Antibody

### Description

<b>Product type</b>	Antibody
<b>Code</b>	BT-MCA3710
<b>Host</b>	Mouse
<b>Isotype</b>	Mouse IgG1
<b>Size</b>	100 $\mu$ L, 50 $\mu$ L
<b>Immunogen</b>	Purified recombinant fragment of human D117L (AA: 1-117) expressed in E. Coli.
<b>Mol wt</b>	N/A
<b>Species reactivity</b>	Others
<b>Clonality</b>	Monoclonal
<b>Recommended application</b>	WB
<b>Concentration</b>	N/A
<b>Full name</b>	N/A
<b>Synonyms</b>	p17

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

### Background


structural protein p17 [ African swine fever virus ]

### Recommended Dilution

WB: 1:500 - 1:2000

Not yet tested in other applications.

### Images



Western blot analysis using D117L mAb against human D117L(AA: 1-117) recombinant protein. lane 1 :(100 ng); lane 2 :(50 ng); lane 3 :(25 ng); lane 4 :(10 ng); lane 5 :(2.5 ng); (Expected MW is 40 kDa)

### Storage

Store at 4°C short term. Aliquot and store at -20°C long term.

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: [save@bt-laboratory.com](mailto:save@bt-laboratory.com) | [www.bt-laboratory.com](http://www.bt-laboratory.com)