

TFAP2B Monoclonal Antibody

Description

Product type	Antibody
Code	BT-MCA3744
Host	Mouse
Isotype	Mouse IgG2b
Size	100µL, 50µL
Immunogen	Purified recombinant fragment of human TFAP2B (AA: 84-193) expressed in E. Coli.
Mol wt	50.5kDa
Species reactivity	Human
Clonality	Monoclonal
Recommended application	WB
Concentration	N/A
Full name	N/A
Synonyms	PDA2;AP-2B;AP2-B

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

N/A


Recommended Dilution

WB: 1:500 - 1:2000


ELISA: 1:10000

Not yet tested in other applications.


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using TFAP2B mAb against human TFAP2B (AA: 84-193) recombinant protein. (Expected MW is 38 kDa)



Western blot analysis using TFAP2B mAb against HEK293 (1) and TFAP2B (AA: 84-193)-hlgGfc transfected HEK293 (2) cell lysate.



Western blot analysis using TFAP2B mouse mAb against SK-N-SH (1) cell lysate.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term.

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com