

## PAX3 Monoclonal Antibody

### Description

<b>Product type</b>	Antibody
<b>Code</b>	BT-MCA3778
<b>Host</b>	Mouse
<b>Isotype</b>	Mouse IgG1
<b>Size</b>	100µL, 50µL
<b>Immunogen</b>	Purified recombinant fragment of human PAX3 (AA: 142-203) expressed in E. Coli.
<b>Mol wt</b>	53kDa
<b>Species reactivity</b>	Human,Mouse
<b>Clonality</b>	Monoclonal
<b>Recommended application</b>	WB
<b>Concentration</b>	N/A
<b>Full name</b>	N/A
<b>Synonyms</b>	WS1;WS3;CDHS;HUP2

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

### Background

N/A


### Recommended Dilution

WB: 1:500 - 1:2000


ELISA: 1:10000

Not yet tested in other applications.


### Images




Western blot analysis using PAX3 mAb against human PAX3 recombinant protein. (Expected MW is 32.6 kDa)



Western blot analysis using PAX3 mAb against HEK293 (1) and PAX3 (AA: 142-203)-hIgGfC transfected HEK293 (2) cell lysate.



Western blot analysis using PAX3 mouse mAb against A549 (1), A431 (2), Jurkat (3), Rat spleen (4) and Mouse liver (5) cell lysate.



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);

### Storage

Store at 4°C short term. Aliquot and store at -20°C long term.

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: [save@bt-laboratory.com](mailto:save@bt-laboratory.com) | [www.bt-laboratory.com](http://www.bt-laboratory.com)