

# Ractopamine Monoclonal Antibody

## Description

Product type Antibody

Code BT-MCA3803

Host Mouse

Isotype Mouse IgG

Size  $100\mu L$ ,  $50\mu L$ 

Immunogen Rac-OVA

Mol wt N/A

Species reactivity Others

Clonality Monoclonal

Recommended application Others

Concentration N/A

 $\label{eq:full name N/A} \textbf{Synonyms} \qquad \qquad \textbf{N/A}$   $\textbf{Synonyms} \qquad \qquad \textbf{N/A}$ 

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

### Background

Ractopamine is a phenolethanolamine

# Recommended Dilution

ELISA: 1:10000

Not yet tested in other applications.

#### Images

No images.

### Storage

Store at 4°C short term. Aliquot and store at -20°C long term.

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com