

GAB2 Monoclonal Antibody

Description

Product type Antibody

Code BT-MCA3941

Host Mouse

 Isotype
 Mouse IgG2a

 Size
 100μL, 50μL

Immunogen Purified recombinant fragment of human GAB2 (AA: 332-476) expressed in E. Coli.

Mol wt 74.5kDa

Species reactivity Human

Clonality Monoclonal

Recommended application IHC,ICC,FCM

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background


This gene is a member of the GRB2-associated binding protein (GAB) gene family. These proteins contain pleckstrin homology (PH) domain, and bind SHP2 tyrosine phosphatase and GRB2 adapter protein. They act as adapters for transmitting various signals in response to stimuli through cytokine and growth factor receptors, and T- and B-cell antigen receptors. The protein encoded by this gene is the principal activator of phosphatidylinositol-3 kinase in response to activation of the high affinity IgE receptor. Two alternatively spliced transcripts encoding different isoforms have been described for this gene.

Recommended Dilution

WB: 1:500 - 1:2000 IHC-p: 1:200 - 1:1000 ICC: 1:200 - 1:1000 FCM: 1:200 - 1:400 ELISA: 1:10000


Not yet tested in other applications.

Images




Western blot analysis using GAB2 mAb against human GAB2 (AA: 332-476) recombinant protein. (Expected MW is 41 kDa)


Western blot analysis using GAB2 mAb against HEK293 (1) and GAB2 (AA: 332-476)-hIgGFc transfected HEK293 (2) cell lysate.




Immunofluorescence analysis of Hela cells using GAB2 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Secondary antibody from Fisher (Cat#: 35503)




Flow cytometric analysis of Hela cells using GAB2 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded cervical cancer tissues using GAB2 mouse mAb with DAB staining.



 $Immun ohistochemical \ analysis \ of paraffin-embedded \ esophageal \ cancer \ tissues \ using \ GAB2 \ mouse$ $mAb \ with \ DAB \ staining.$



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);

Storage

Store at 4°C short term. Aliquot and store at -20°C long term.

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com