

MCAM Monoclonal Antibody

Description

Product type	Antibody
Code	BT-MCA3996
Host	Mouse
Isotype	Mouse IgG1
Size	100µL, 50µL
Immunogen	Purified recombinant fragment of human MCAM (AA: 84-189) expressed in E. Coli.
Mol wt	71.6kDa
Species reactivity	Human
Clonality	Monoclonal
Recommended application	WB,IHC
Concentration	N/A
Full name	N/A
Synonyms	CD146;MUC18

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

N/A

Recommended Dilution


WB: 1:500 - 1:2000

IHC-p: 1:200 - 1:1000


ELISA: 1:10000

Not yet tested in other applications.


Images




Western blot analysis using MCAM mAb against human MCAM recombinant protein. (Expected MW is 37.7 kDa)




Western blot analysis using MCAM mAb against HEK293 (1) and MCAM (AA: 84-189)-hIgGfC transfected HEK293 (2) cell lysate.




Western blot analysis using MCAM mouse mAb against Hela cell lysate.



Immunohistochemical analysis of paraffin-embedded liver cancer tissues using MCAM mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded kidney tissues using MCAM mouse mAb with DAB staining.



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);

Storage

Store at 4°C short term. Aliquot and store at -20°C long term.

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com