

ELK1 Monoclonal Antibody

Description

Product type	Antibody
Code	BT-MCA4085
Host	Mouse
Isotype	Mouse IgG1
Size	100μL, 50μL
Immunogen	Purified recombinant fragment of ELK1 expressed in E. Coli.
Mol wt	N/A
Species reactivity	Human
Clonality	Monoclonal
Recommended application	IHC
Concentration	N/A
Full name	N/A
Synonyms	ELK1

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.


Background

N/A

Recommended Dilution

WB: 1:500 - 1:2000 IHC-p: 1:200 - 1:1000 ELISA: 1:10000 Not yet tested in other applications.

Images



Western blot analysis using ELK1 mouse mAb against truncated ELK1 recombinant protein (1) and K562 cell lysate (2).

Immunohistochemical analysis of paraffin-embedded human brain tumor tissue, showing nuclear and cytoplasmic localization using ELK1 mouse mAb with DAB staining.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term.

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com