

ULK2 Monoclonal Antibody

Description

Product type	Antibody
Code	BT-MCA4207
Host	Mouse
Isotype	Mouse IgG1
Size	100μL, 50μL
Immunogen	Purified recombinant fragment of human ULK2 (AA: 1-155) expressed in E. Coli.
Mol wt	112.7kDa
Species reactivity	Human,Mouse,Rat
Clonality	Monoclonal
Recommended application	WB,FCM
Concentration	N/A
Full name	N/A
Synonyms	ATG1B;Unc51.2

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.


Background


N/A

Recommended Dilution

WB: 1:500 - 1:2000 FCM: 1:200 - 1:400 ELISA: 1:10000 Not yet tested in other applications.

Images






Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)

Western blot analysis using ULK2 mAb against human ULK2 (AA: 1-155) recombinant protein. (Expected MW is 43.4 kDa)

Western blot analysis using ULK2 mAb against HEK293 (1) and ULK2 (AA: 1-155)-hIgGFc transfected HEK293 (2) cell lysate.



Western blot analysis using ULK2 mouse mAb against NIH/3T3 (1), HepG2 (2), SK-Hep-1 (3), SK-OV-3 (4), C6 (5), PC-12 (6), and MCF-7 (7) cell lysate.

Flow cytometric analysis of Hela cells using ULK2 mouse mAb (green) and negative control (red).

Storage

Store at 4°C short term. Aliquot and store at -20°C long term.

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com