

OX40 Monoclonal Antibody

Description

Product type	Antibody
Code	BT-MCA4291
Host	Mouse
Isotype	Mouse IgG2a
Size	100μL, 50μL
Immunogen	Purified recombinant fragment of human OX40 (AA: extra 29-214) expressed in HEK293 cells.
Mol wt	29.3kDa
Species reactivity	Human
Clonality	Monoclonal
Recommended application	Others
Concentration	N/A
Full name	N/A
Synonyms	TNFRSF4;ACT35;CD134;IMD16;TXGP1L

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

N/A


Recommended Dilution

WB: 1:500 - 1:2000


ELISA: 1:10000

Not yet tested in other applications.

Images



Black line: Control Antigen (100 ng); Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)



Western blot analysis using OX40 mAb against human OX40 (AA: extra 29-214) recombinant protein. (Expected MW is 50 kDa)

Storage

Store at 4°C short term. Aliquot and store at -20°C long term.