

# HAVCR1 Monoclonal Antibody

## Description

<b>Product type</b>	Antibody
<b>Code</b>	BT-MCA4631
<b>Host</b>	Mouse
<b>Isotype</b>	Mouse IgG1
<b>Size</b>	100µL, 50µL
<b>Immunogen</b>	Purified recombinant fragment of human HAVCR1 (AA: 70-290) expressed in E. Coli.
<b>Mol wt</b>	38.7kDa
<b>Species reactivity</b>	Human,Mouse
<b>Clonality</b>	Monoclonal
<b>Recommended application</b>	WB
<b>Concentration</b>	N/A
<b>Full name</b>	N/A
<b>Synonyms</b>	TIM;KIM1;TIM1;CD365;HAVCR;KIM-1;TIM-1;TIMD1;TIMD-1;HAVCR-1

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

## Background

N/A


## Recommended Dilution

WB: 1:500 - 1:2000


ELISA: 1:10000

Not yet tested in other applications.


## Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using HAVCR1 mAb against human HAVCR1 (AA: 70-290) recombinant protein. (Expected MW is 37 kDa)



Western blot analysis using HAVCR1 mAb against HEK293 (1) and HAVCR1 (AA: 70-290)-hlgGfC transfected HEK293 (2) cell lysate.



Western blot analysis using HAVCR1 mouse mAb against NIH/3T3 (1), HEK293 (2), HeLa (3), Raw264.7 (4), Jurkat (5), and PC-2 (6) cell lysate.

### Storage

Store at 4°C short term. Aliquot and store at -20°C long term.

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: [save@bt-laboratory.com](mailto:save@bt-laboratory.com) | [www.bt-laboratory.com](http://www.bt-laboratory.com)