

PDGFR Monoclonal Antibody

Description

Product type	Antibody
Code	BT-MCA4657
Host	Mouse
Isotype	Mouse IgG1
Size	100 μ L, 50 μ L
Immunogen	Purified recombinant fragment of human PDGFR β expressed in E. Coli.
Mol wt	190kDa
Species reactivity	Human,Mouse
Clonality	Monoclonal
Recommended application	WB
Concentration	N/A
Full name	N/A
Synonyms	JTK12;PDGFR;CD140B;PDGFR1;PDGFRB

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Platelet-derived growth factor receptor, beta polypeptide, also known as PDGFRB, PDGFR


Recommended Dilution

WB: 1:500 - 1:2000

ELISA: 1:10000

Not yet tested in other applications.

Images



Western blot analysis using PDGFRb mouse mAb against NIH/3T3 cell lysate (1).

Storage

Store at 4°C short term. Aliquot and store at -20°C long term.

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com