

HSPA9 Monoclonal Antibody

Description

Product type	Antibody
Code	BT-MCA4694
Host	Mouse
Isotype	Mouse IgG1
Size	100μL, 50μL
Immunogen	Purified recombinant fragment of human HSPA9 (AA: 480-679) expressed in mammalian.
Mol wt	74KDa
Species reactivity	Human,Monkey,Rat
Clonality	Monoclonal
Recommended application	WB,IHC,ICC,FCM
Concentration	N/A
Full name	N/A
Synonyms	CSA;MOT;MOT2;SAAN;CRP40;EVPLS;GRP75;PBP74;GRP-75;HSPA9B;SIDBA4;MTHSP75;HEL-S-124m
	124m

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.


Background

This gene encodes a member of the heat shock protein 70 gene family. The encoded protein is primarily localized to the mitochondria but is also found in the endoplasmic reticulum, plasma membrane and cytoplasmic vesicles. This protein is a heat-shock cognate protein. This protein plays a role in cell proliferation, stress response and maintenance of the mitochondria. A pseudogene of this gene is found on chromosome 2.

Recommended Dilution


WB: 1:500 - 1:2000 IHC-p: 1:200 - 1:1000 ICC: 1:200 - 1:1000 FCM: 1:200 - 1:400 ELISA: 1:10000 Not yet tested in other applications.

Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)

Western blot analysis using HSPA9 mAb against human HSPA9 (AA: 480-679) recombinant protein. (Expected MW is 25.2 kDa)




Immunofluorescence analysis of Hela cells using HSPA9 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)

Western blot analysis using HSPA9 mouse mAb against CSO-7 (1),C6 (2), PC-12 (3), PANC-1(4),A549 (5),MCF-7 (6),K562 (7),Hela (8),A431 (9),HepG2 (10)and Jurkat (11) cell lysate.

Flow cytometric analysis of Hela cells using HSPA9 mouse mAb (green) and negative control (red).

Flow cytometric analysis of HepG2 cells using HSPA9 mouse mAb (green) and negative control (red).

Flow cytometric analysis of Jurkat cells using HSPA9 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded cervical carcinoma tissues using HSPA9 mouse mAb with DAB staining.

Immunohistochemical analysis of paraffin-embedded ovarian cancer tissues using HSPA9 mouse mAb with DAB staining.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term.

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com