

CBX8 Monoclonal Antibody

Description

Product type	Antibody
Code	BT-MCA4809
Host	Mouse
Isotype	Mouse IgG1
Size	100µL, 50µL
Immunogen	Purified recombinant fragment of human CBX8 (AA: 17-222) expressed in E. Coli.
Mol wt	43.4kDa
Species reactivity	Human
Clonality	Monoclonal
Recommended application	WB,FCM
Concentration	N/A
Full name	N/A
Synonyms	PC3;RC1

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

N/A

Recommended Dilution


WB: 1:500 - 1:2000

FCM: 1:200 - 1:400


ELISA: 1:10000

Not yet tested in other applications.


Images




Western blot analysis using CBX8 mAb against human CBX8 recombinant protein. (Expected MW is 49.5 kDa)



Western blot analysis using CBX8 mAb against HEK293 (1) and CBX8 (AA: 17-222)-hIgGfc transfected HEK293 (2) cell lysate.



Flow cytometric analysis of HeLa cells using CBX8 mouse mAb (green) and negative control (red).



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);

Storage

Store at 4°C short term. Aliquot and store at -20°C long term.

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com