

PKHD1 Monoclonal Antibody

Description

Product type	Antibody
Code	BT-MCA4815
Host	Rat
Isotype	Rat IgG1
Size	100µL, 50µL
Immunogen	Purified recombinant fragment of mouse PKHD1 (AA: 3878-4060) expressed in E. Coli.
Mol wt	445kDa
Species reactivity	Human
Clonality	Monoclonal
Recommended application	IHC,ICC,FCM
Concentration	N/A
Full name	N/A
Synonyms	FPC;Tigm1;AI118496;AI182499

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

N/A

Recommended Dilution

IHC-p: 1:200 - 1:1000


ICC: 1:200 - 1:1000

FCM: 1:200 - 1:400


ELISA: 1:10000

Not yet tested in other applications.


Images




Western blot analysis using PKHD1 mAb against mouse PKHD1(AA: 3878-4060) recombinant protein. (Expected MW is 23 kDa)




Immunofluorescence analysis of A431 cells using PKHD1 Rat mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)



Flow cytometric analysis of HeLa cells using PKHD1 Rat mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded esophageal cancer tissues using PKHD1 Rat mAb with DAB staining.



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);

Storage

Store at 4°C short term. Aliquot and store at -20°C long term.

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com